



DROSHA Antibody, HRP conjugated

Product Code	CSB-PA889096LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRR4
Immunogen	Recombinant Human Ribonuclease 3 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease 3 (EC 3.1.26.3) (Protein Drosha) (Ribonuclease III) (RNase III) (p241), DROSHA, RN3 RNASE3L RNASEN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DROSHA