



DLG4 Antibody, HRP conjugated

Product Code	CSB-PA006938LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78352
Immunogen	Recombinant Human Disks large homolog 4 protein (551-724AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disks large homolog 4 (Postsynaptic density protein 95) (PSD-95) (Synapse-associated protein 90) (SAP-90) (SAP90), DLG4, PSD95
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DLG4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.