



GRK7 Antibody? Biotin conjugated

Product Code	CSB-PA009929LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTQ7
Immunogen	Recombinant Human Rhodopsin kinase protein (342-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhodopsin kinase (EC 2.7.11.14) (G protein-coupled receptor kinase 7) (G protein-coupled receptor kinase GRK7), GRK7, GPRK7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GRK7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.