



KCNK15 Antibody, HRP conjugated

Product Code	CSB-PA867147LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H427
Immunogen	Recombinant Human Potassium channel subfamily K member 15 protein (244-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel subfamily K member 15 (Acid-sensitive potassium channel protein TASK-5) (TWIK-related acid-sensitive K(+) channel 5) (Two pore potassium channel KT3.3) (Two pore K(+) channel KT3.3), KCNK15, TASK5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KCNK15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.