



# Ace2 Antibody, HRP conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA688690OB01RA  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q5EGZ1  |
| <b>Immunogen</b>           | Peptide sequence from Rat Ace2 protein (105-118AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Rat   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | HRP   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| <b>Purification Method</b> | >95%, Antigen Affinity purified & Affinity purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Angiotensin-converting enzyme 2 (EC 3.4.17.23) (ACE-related carboxypeptidase) [Cleaved into: Processed angiotensin-converting enzyme 2], Ace2 |
| <b>Immunogen Species</b>   | Rattus norvegicus (Rat)   |
| <b>Research Area</b>       | Signal Transduction   |
| <b>Target Names</b>        | Ace2  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |