



Ace2 Antibody, FITC conjugated

Product Code	CSB-PA688690OC01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5EGZ1
Immunogen	Peptide sequence from Rat Ace2 protein (105-118AA)
Raised In	Rabbit
Species Reactivity	Rat
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiotensin-converting enzyme 2 (EC 3.4.17.23) (ACE-related carboxypeptidase) [Cleaved into: Processed angiotensin-converting enzyme 2], Ace2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Ace2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.