



AGO1 Antibody, Biotin conjugated

Product Code	CSB-PA892141LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL18
Immunogen	Recombinant Human Protein argonaute-1 protein (251-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein argonaute-1 (Argonaute1) (hAgo1) (Argonaute RISC catalytic component 1) (Eukaryotic translation initiation factor 2C 1) (eIF-2C 1) (eIF2C 1) (Putative RNA-binding protein Q99), AGO1, EIF2C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AGO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.