




# Apoa1 Antibody

<b>Product Code</b>	CSB-PA001913LA01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04639
<b>Immunogen</b>	Recombinant Rat ApoA1 protein (19-259AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Antigen Affinity purified & Affinity purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Apolipoprotein A-I (Apo-AI) (ApoA-I) (Apolipoprotein A1) [Cleaved into: Proapolipoprotein A-I (ProapoA-I)], ApoA1
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	Apoa1

## Image



### Western Blot

Positive WB detected in: Rat Heart tissue, Rat Liver tissue, Rat Spleen tissue, Rat Lung tissue, Rat Kidney tissue, Rat Brain tissue, Rat Stomach tissue

All lanes: ApoA1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.