



APOO Antibody, FITC conjugated

Product Code	CSB-PA001948LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUR5
Immunogen	Recombinant Human MICOS complex subunit MIC26 protein (26-198AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MICOS complex subunit MIC26 (Apolipoprotein O) (MICOS complex subunit MIC23) (Protein FAM121B), APOO, FAM121B MIC23 MIC26
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	APOO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.