



ATF6 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA002277LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P18850 |
| Immunogen | Recombinant Human Cyclic AMP-dependent transcription factor ATF-6 alpha protein (1-194AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Antigen Affinity purified & Affinity purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cyclic AMP-dependent transcription factor ATF-6 alpha (cAMP-dependent transcription factor ATF-6 alpha) (Activating transcription factor 6 alpha) (ATF6-alpha) [Cleaved into: Processed cyclic AMP-dependent transcription factor ATF-6 alpha], ATF6 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ATF6 |