



ATF6 Antibody, Biotin conjugated

Product Code	CSB-PA002277LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18850
Immunogen	Recombinant Human Cyclic AMP-dependent transcription factor ATF-6 alpha protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclic AMP-dependent transcription factor ATF-6 alpha (cAMP-dependent transcription factor ATF-6 alpha) (Activating transcription factor 6 alpha) (ATF6-alpha) [Cleaved into: Processed cyclic AMP-dependent transcription factor ATF-6 alpha], ATF6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.