



BHLHE23 Antibody, Biotin conjugated

Product Code	CSB-PA850867LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDY6
Immunogen	Recombinant Human Class E basic helix-loop-helix protein 23 protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Class E basic helix-loop-helix protein 23 (bHLHe23) (Class B basic helix-loop-helix protein 4) (bHLHb4), BHLHE23, BHLHB4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BHLHE23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.