



CALR Antibody, Biotin conjugated

Product Code	CSB-PA004458OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27797
Immunogen	Peptide sequence from Human Calreticulin protein (185-196AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calreticulin (CRP55) (Calregulin) (Endoplasmic reticulum resident protein 60) (ERp60) (HACBP) (grp60), CALR, CRTC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CALR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.