



CAMK2A Antibody, Biotin conjugated

Product Code	CSB-PA892176HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQM7
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase type II subunit alpha protein (321-478AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase type II subunit alpha (CaM kinase II subunit alpha) (CaMK-II subunit alpha) (EC 2.7.11.17), CAMK2A, CAMKA KIAA0968
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAMK2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.