



CDK6 Antibody, Biotin conjugated

Product Code	CSB-PA005074LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00534
Immunogen	Recombinant Human Cyclin-dependent kinase 6 protein (102-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 6 (EC 2.7.11.22) (Cell division protein kinase 6) (Serine/threonine-protein kinase PLSTIRE), CDK6, CDKN6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CDK6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.