



CNR1 Antibody, Biotin conjugated

Product Code	CSB-PA005678LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21554
Immunogen	Recombinant Human Cannabinoid receptor 1 protein (1-116AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cannabinoid receptor 1 (CB-R) (CB1) (CANN6), CNR1, CNR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CNR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.