



CNTNAP3B Antibody, Biotin conjugated

Product Code	CSB-PA822273LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NU0
Immunogen	Recombinant Human Contactin-associated protein-like 3B protein (26-256AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Contactin-associated protein-like 3B (Cell recognition molecule Caspr3b), CNTNAP3B, CASPR3B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CNTNAP3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.