



COL18A1 Antibody, Biotin conjugated

Product Code	CSB-PA005725LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39060
Immunogen	Recombinant Human Collagen alpha-1(XVIII) chain protein (1600-1754AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-1(XVIII) chain [Cleaved into: Endostatin, Non-collagenous domain 1 (NC1)], COL18A1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	COL18A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.