



FAM24B Antibody, Biotin conjugated

Product Code	CSB-PA843275LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5W8
Immunogen	Recombinant Human Protein FAM24B protein (22-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM24B, FAM24B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAM24B