



FBXO36 Antibody, FITC conjugated

Product Code	CSB-PA854107LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEA4
Immunogen	Recombinant Human F-box only protein 36 protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 36, FBXO36, FBX36
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FBXO36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.