




FST Antibody

Product Code	CSB-PA009024LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19883
Immunogen	Recombinant Human Follistatin protein (30-344AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Follistatin (FS) (Activin-binding protein), FST
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FST

Image



Western Blot

Positive WB detected in: U251 whole cell lysate

All lanes: FST antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39, 35 kDa

Observed band size: 39 kDa