




GCG Antibody

Product Code	CSB-PA009315LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01275
Immunogen	Recombinant Human Glucagon protein (21-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucagon [Cleaved into: Glicentin, Glicentin-related polypeptide (GRPP), Oxyntomodulin (OXM) (OXY), Glucagon, Glucagon-like peptide 1 (GLP-1) (Incretin hormone), Glucagon-like peptide 1(7-37) (GLP-1(7-37)), Glucagon-like peptide 1(7-36) (GLP-1(7-36)), Glucagon-like peptide 2 (GLP-2)], GCG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GCG

Image



IHC image of CSB-PA009315LA01HU diluted at 1:300 and staining in paraffin-embedded human pancreatic tissue performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a Goat anti-rabbit polymer IgG labeled by HRP and visualized using 0.05% DAB.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.