



GDF11 Antibody, Biotin conjugated

Product Code	CSB-PA009344LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95390
Immunogen	Recombinant Human Growth/differentiation factor 11 protein (299-407AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth/differentiation factor 11 (GDF-11) (Bone morphogenetic protein 11) (BMP-11), GDF11, BMP11
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GDF11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.