




GLB1L3 Antibody

Product Code	CSB-PA850859LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCI6
Immunogen	Recombinant Human Beta-galactosidase-1-like protein 3 protein (1-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-galactosidase-1-like protein 3 (EC 3.2.1.-), GLB1L3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GLB1L3

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate, HeLa whole cell lysate

All lanes: GLB1L3 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75, 36 kDa

Observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.