



GNLY Antibody, FITC conjugated

Product Code	CSB-PA009627LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22749
Immunogen	Recombinant Human Granulysin protein (23-145AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granulysin (Lymphokine LAG-2) (Protein NKG5) (T-cell activation protein 519), GNLY, LAG2 NKG5 TLA519
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GNLY