



GPLY Antibody, Biotin conjugated

Product Code	CSB-PA009627LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22749
Immunogen	Recombinant Human Granulysin protein (23-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granulysin (Lymphokine LAG-2) (Protein NKG5) (T-cell activation protein 519), GPLY, LAG2 NKG5 TLA519
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPLY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.