




Hmgb1 Antibody

Product Code	CSB-PA010553NA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63159
Immunogen	Peptide sequence from Rat Hmgb1 protein (2-18AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	High mobility group protein B1 (Amphoterin) (Heparin-binding protein p30) (High mobility group protein 1) (HMG-1), Hmgb1, Hmg-1 Hmg1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Hmgb1

Image



Western Blot

Positive WB detected in: Rat Brain tissue

All lanes: Hmgb1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 25 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.