



Hmgb1 Antibody, Biotin conjugated

Product Code	CSB-PA010553ND01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63159
Immunogen	Peptide sequence from Rat Hmgb1 protein (2-18AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	High mobility group protein B1 (Amphoterin) (Heparin-binding protein p30) (High mobility group protein 1) (HMG-1), Hmgb1, Hmg-1 Hmg1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Hmgb1