



IGHG3 Antibody, HRP conjugated

Product Code	CSB-PA011158OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01860
Immunogen	Peptide sequence from Human Immunoglobulin heavy constant gamma 3 protein (312-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin heavy constant gamma 3 (HDC) (Heavy chain disease protein) (Ig gamma-3 chain C region), IGHG3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IGHG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.