



Lbp Antibody, Biotin conjugated

Product Code	CSB-PA714125OD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63313
Immunogen	Peptide sequence from Rat Lipopolysaccharide-binding protein (250-264AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipopolysaccharide-binding protein (LBP), Lbp
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Lbp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.