



Lcn2 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA012816LB01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P30152 |
| Immunogen | Recombinant Rat Lcn2 protein (1-198AA) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Neutrophil gelatinase-associated lipocalin (NGAL) (Alpha-2-microglobulin-related protein) (Alpha-2U globulin-related protein) (Lipocalin-2) (Siderocalin LCN2) (p25), Lcn2 |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Cardiovascular |
| Target Names | Lcn2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |