



Lcn2 Antibody, Biotin conjugated

Product Code	CSB-PA012816LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30152
Immunogen	Recombinant Rat Lcn2 protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil gelatinase-associated lipocalin (NGAL) (Alpha-2-microglobulin-related protein) (Alpha-2U globulin-related protein) (Lipocalin-2) (Siderocalin LCN2) (p25), Lcn2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Lcn2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.