



Mmp2 Antibody, HRP conjugated

Product Code	CSB-PA014666LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33434
Immunogen	Recombinant Mouse 72 kDa type IV collagenase protein (412-662AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	72 kDa type IV collagenase (EC 3.4.24.24) (72 kDa gelatinase) (Gelatinase A) (Matrix metalloproteinase-2) (MMP-2) [Cleaved into: PEX], Mmp2
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Mmp2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.