




# MROH8 Antibody

<b>Product Code</b>	CSB-PA872499LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9H579
<b>Immunogen</b>	Recombinant Human Protein MROH8 protein (1-275AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein MROH8 (Maestro heat-like repeat-containing protein family member 8), MROH8, C20orf131 C20orf132
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	MROH8

## Image



### Western Blot

Positive WB detected in: U87 whole cell lysate  
 All lanes: MROH8 antibody at 1:2000  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/50000 dilution  
 Predicted band size: 55, 100 kDa  
 Observed band size: 55 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.