



MROH8 Antibody, HRP conjugated

Product Code	CSB-PA872499LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H579
Immunogen	Recombinant Human Protein MROH8 protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein MROH8 (Maestro heat-like repeat-containing protein family member 8), MROH8, C20orf131 C20orf132
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MROH8