



MROH8 Antibody, Biotin conjugated

Product Code	CSB-PA872499LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H579
Immunogen	Recombinant Human Protein MROH8 protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein MROH8 (Maestro heat-like repeat-containing protein family member 8), MROH8, C20orf131 C20orf132
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MROH8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.