




Nampt Antibody

Product Code	CSB-PA015422LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q80Z29
Immunogen	Recombinant Rat Nampt protein (231-491AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotinamide phosphoribosyltransferase (NAMPTase) (Nampt) (EC 2.4.2.12) (Pre-B-cell colony-enhancing factor 1 homolog) (PBEF) (Visfatin), Nampt, Pbef1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Nampt

Image



Western Blot

Positive WB detected in: Rat Heart tissue, Rat Liver tissue, Rat Spleen tissue, Rat Kidney tissue, Rat Brain tissue, HeLa whole cell lysate
All lanes: Nampt antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 56 kDa
Observed band size: 56 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.