



# Nampt Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA015422LB01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q80Z29
<b>Immunogen</b>	Recombinant Rat Nampt protein (231-491AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Antigen Affinity purified & Affinity purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Nicotinamide phosphoribosyltransferase (NAMPRase) (Nampt) (EC 2.4.2.12) (Pre-B-cell colony-enhancing factor 1 homolog) (PBEF) (Visfatin), Nampt, Pbef1
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	Nampt
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.