



Nampt Antibody, Biotin conjugated

Product Code	CSB-PA015422LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q80Z29
Immunogen	Recombinant Rat Nampt protein (231-491AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotinamide phosphoribosyltransferase (NAMPRtase) (Nampt) (EC 2.4.2.12) (Pre-B-cell colony-enhancing factor 1 homolog) (PBEF) (Visfatin), Nampt, Pbef1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Nampt
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.