



# NFATC2 Antibody, HRP conjugated

Product Code	CSB-PA615689LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13469
Immunogen	Recombinant Human Nuclear factor of activated T-cells, cytoplasmic 2 protein (210-440AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor of activated T-cells, cytoplasmic 2 (NF-ATc2) (NFATc2) (NFAT pre-existing subunit) (NF-ATp) (T-cell transcription factor NFAT1), NFATC2, NFAT1 NFATP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NFATC2