



NFATC2 Antibody, Biotin conjugated

Product Code	CSB-PA615689LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13469
Immunogen	Recombinant Human Nuclear factor of activated T-cells, cytoplasmic 2 protein (210-440AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor of activated T-cells, cytoplasmic 2 (NF-ATc2) (NFATc2) (NFAT pre-existing subunit) (NF-ATp) (T-cell transcription factor NFAT1), NFATC2, NFAT1 NFATP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NFATC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.