



NFE2L2 Antibody, FITC conjugated

Product Code	CSB-PA614961LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16236
Immunogen	Recombinant Human Nuclear factor erythroid 2-related factor 2 protein (330-448AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor erythroid 2-related factor 2 (NF-E2-related factor 2) (NFE2-related factor 2) (HEBP1) (Nuclear factor, erythroid derived 2, like 2), NFE2L2, NRF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NFE2L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.