



NFKBIA Antibody, HRP conjugated

Product Code	CSB-PA015761LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25963
Immunogen	Recombinant Human NF-kappa-B inhibitor alpha protein (27-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	NF-kappa-B inhibitor alpha (I-kappa-B-alpha) (Ikb-alpha) (IkappaBalpha) (Major histocompatibility complex enhancer-binding protein MAD3), NFKBIA, IKBA MAD3 NFKBI
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NFKBIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.