




Nucb2 Antibody

Product Code	CSB-PA884912LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JI85
Immunogen	Recombinant Rat Nucb2 protein (194-420AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleobindin-2 (DNA-binding protein NEFA) (Prepronesfatin) [Cleaved into: Nesfatin-1], Nucb2, Nefa
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Nucb2

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate

All lanes: Nucb2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51 kDa

Observed band size: 51 kDa