



OR10A2 Antibody, Biotin conjugated

Product Code	CSB-PA884445ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H208
Immunogen	Peptide sequence from Human Olfactory receptor 10A2 protein (247-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 10A2 (HP4) (Olfactory receptor OR11-86), OR10A2, OR10A2P
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR10A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.