




OR10P1 Antibody

Product Code	CSB-PA818792OA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGE3
Immunogen	Peptide sequence from Human Olfactory receptor 10P1 protein (260-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 10P1 (Olfactory receptor 10P2) (Olfactory receptor 10P3) (Olfactory receptor OR12-7), OR10P1, OR10P1P OR10P2P OR10P3P
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR10P1

Image



Western Blot

Positive WB detected in: HL60 whole cell lysate

All lanes: OR10P1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa