



OR1K1 Antibody, Biotin conjugated

Product Code	CSB-PA847729ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGR3
Immunogen	Peptide sequence from Human Olfactory receptor 1K1 protein (261-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 1K1, OR1K1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR1K1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.