



OR2A7 Antibody, Biotin conjugated

Product Code	CSB-PA857015OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96R45
Immunogen	Peptide sequence from Human Olfactory receptor 2A7 protein (258-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 2A7 (Olfactory receptor OR7-18), OR2A7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR2A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.