



OR2M2 Antibody, HRP conjugated

Product Code	CSB-PA850420OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96R28
Immunogen	Peptide sequence from Human Olfactory receptor 2M2 protein (259-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 2M2 (OST423), OR2M2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR2M2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.