



OR2S2 Antibody, Biotin conjugated

Product Code	CSB-PA873623OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQN1
Immunogen	Peptide sequence from Human Olfactory receptor 2S2 protein (262-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 2S2 (Olfactory receptor OR9-3), OR2S2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OR2S2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.